


Capture Compound Mass Spectrometry®: on- and off-target deconvolution

T. Ladduwaheetty, S. Almond*, E. Bush*, A. Cridland*, S. Dowler*, G. Hardman*, D. Kenny*, I. Linney*, N. Macabuag*, D. Mitchell*, P. Mitchell*
 Charles River Laboratories, Chesterford Research Park, Saffron Walden, UK
 Tammy.Ladduwaheetty@crl.com
 * CCMS team at Charles River; Sarah.Almond@crl.com



1 CCMS INTRODUCTION

- Capture Compound Mass Spectrometry (CCMS) is an unbiased, proteome-wide approach for the identification of specific on- and off-target binding protein targets for small molecules of interest.
- Capture Compounds® (CCs) are unique bi-functional probes designed to interrogate proteins in their native environments. The distinct molecular architecture of the Capture Compounds® enables a three-stage process of binding.



2 CCMS is a chemoproteomic tool for profiling protein interactions of small molecules


- Unbiased screening of the entire proteome compared with traditional selectivity screening



CCMS can support projects at a variety of stages throughout drug discovery and development.


3 THE CCMS EXPERIMENT

- A panel of CCs are synthesized with the selectivity function in a position compliant with on-target SAR and in alternative orientations to allow generation of a full interaction profile
- CCs are incubated with a biological sample whereby binding of the selectivity function to interacting proteins occurs before covalent capture via photo-irradiation
- Competition experiments include incubation of the sample with both the CC and an excess of the free ligand to allow specifically interacting proteins to be determined
- Sample analysis by high-resolution LC-MS/MS reveals specific binding partners



4 CASE STUDY: CCMS IDENTIFIES ON TARGET MECHANISM FROM PHENOTYPIC SCREEN

- Client showed that two structurally distinct series displaying the same phenotype
 - Structural diversity between the two series suggested alternative MoA possible
- A precursor molecule from each series synthesized
 - Each retained activity in the phenotypic assay
- Two Capture Compounds® were synthesized from each series to attempt target deconvolution
- The Capture Compounds® were used with the entire unfractionated lysate and with a membrane enriched fraction



5 CASE STUDY: CCMS IDENTIFIES OFF-TARGET MECHANISM OF TOLCAPONE TOXICITY

- CCMS was used to determine on- and off-target interactions of the catechol-O-methyl transferase (COMT) inhibitor tolcapone in a human liver cancer cell line (HepG2)
- Aim of study: determine hepatotoxicity causing off-target for tolcapone and compare with entacapone which has no toxicity
- Differential binding profiles were seen against on- and off target proteins. These binding profiles correlated with cellular hepatotoxicity

